

# BRIVIEW

## SINGLE USE CHOLEDOCHOSCOPE

- ▶ Tapered-Shape distal end aids insertion
- ▶ Larger working channel available for added flexibility
- ▶ HD image quality on par with conventional scopes

Universal Processor - compatible with range of disposable scopes



Available in 2 sizes for Surgical & Laparoscopic Applications

MODEL	CHV-100 T-E	CHV-US100-H
Field of View	110 °	110 °
Depth of View	5-50 mm	5-50 mm
Deflection up/ Down	180°/130°	120°/120°
Insertion Diameter	5.3 mm	3.0 mm
Distal End Diameter	Ø 5.3 mm	Ø 3.0 mm
Working Length	38.5 cm	67 cm
Working Channel	2.6 mm / 7.5 Fr	1.2 mm / 3.6 Fr

# Benefits of Briview

## Single-Use Choledochoscope

1

**Safety** - A new sterile scope for every patient, reducing the risk of infection

2

**Consistency** - Predicable clinical performance every time - no issue with wear and tear

3

**24 x 7 Availability** - No need to wait for reprocessing of reusable scopes

4

**Cost Transparency** - No Service or Leasing Agreements

5

**Maintenance Free** - No degradation, No reprocessing and No repairs

6

**Compatibility** - works with most commonly used diagnostic and therapeutic tools