


# Balloon In-Deflation Device

Patent Selok® Self-Locking Mechanism with Maximum Pressure of 30atm



Product Code	Product Description	Unit/Box
DID30S	BALLOON IN-DEFLATION DEVICE	1

- Maximum pressure of 30atm
- 20ml clear syringe barrel with smooth moving plunger
- Selok® technology patent - self-locking mechanism for easy locking and Gradual or rapid pressure release
- Ergonomic handle
- Precision manometer
- Optional 3-way stopcock with "ON" or "OFF" handle

DEMAX Balloon 2022 Rev 1

## Medilex Medical

**ROI:** Unit 5 Bromley Business Park, Kilpedder,  
Co. Wicklow A63 X924  
**NI:** Innis House, Hollywood,  
Co. Down, BT18 9HF

[www.medilexmedical.com](http://www.medilexmedical.com)



**KANDO**  
MEDICAL  
a KanDo Medical Company

## Customer Service

**ROI:** +353 1 201 4620  
**NI:** +44 28 9099 5160  
**Email:** [sales@medilexmedical.com](mailto:sales@medilexmedical.com)