

Hemostasis valve (Y connector)



Product Code	Product Description	Unit/Box
DPY02B	HEMOSTASIS VALVE (Y CONNECTOR)	40

- Unique push-click design
 - Minimizes blood loss effectively
 - One-handed operation
 - Complex PCI procedures
- Ergonomic knob design
 - Screw thread surface offers better control
- Large internal diameter facilitates insertion of instruments (2.7±0.1)
- Max. Pressure: 58 psi (Internal standard 180 psi)
- Extension tube optional



DEMAX Y connector 2022 Rev 1

Medilex Medical

ROI: Unit 5 Bromley Business Park, Kilpedder,
Co. Wicklow A63 X924
NI: Innis House, Hollywood,
Co. Down, BT18 9HF

www.medilexmedical.com



KANDO
MEDICAL
a KanDo Medical Company

Customer Service

ROI: +353 1 201 4620
NI: +44 28 9099 5160
Email: sales@medilexmedical.com